

Maxpar® Phospho-Protein Staining

WARNING Before handling any chemicals, refer to the Safety Data Sheet (SDS) provided by the manufacturer, and observe all relevant precautions.

Reagents and Materials

Included in Panel Kit

- Maxpar® Metal Conjugated Antibodies
- Maxpar® Cell Staining Buffer (Cat# 201068)
- Maxpar® Fix I Buffer (5X) (Cat# 201065)
- Maxpar® Fix and Perm Buffer (Cat# 201067)
- Cell-ID™ Intercalator-Ir (Cat# 201192A [125 µM] or 201192B [500 µM])
- Maxpar® Water (Cat# 201069)

Optional

- Cell-ID™ Cisplatin (Cat# 201064)
- Serum-Free and Complete Media
- Fc-receptor Blocking Solution

Other Required Materials and Equipment

- Methanol (Fisher Scientific Cat# BP1105-4)
- Polystyrene or Polypropylene Round-bottom Tubes, 5 mL capacity, 12 x 75 mm
- Polystyrene or Polypropylene Round-bottom Tubes with Cell-Strainer Cap, 5 mL capacity, 12 x 75 mm
- 1.5 ml microfuge tubes
- Pipet tips with aerosol barrier
- Centrifuge capable of holding 5 ml tubes
- Vacuum aspirator

- Vortexer

Important Notes Before Starting

- This protocol should be followed for staining activation-induced phosphorylated antigens. For staining secreted proteins, including cytokines, and other intracellular antigens, please use the Maxpar Cytoplasmic/Secreted Antigen Staining Protocol.
- For cell centrifugation steps, centrifugation should be performed for 5 minutes at 300g before cell fixation, and for 5 minutes at 800g after cell fixation. The increased centrifugation speed after cell fixation will result in greater cell recovery.
- Fluidigm antibodies are pre-titrated and we recommend staining with 1 μ l of each antibody for 3 million cells in a 100 μ l staining volume; however antibodies can be titrated for individual experiments.
- An optional Fc-blocking step is recommended in the following protocol to prevent binding of Maxpar Metal Conjugated Antibodies to Fc receptors, which will result in high non-specific background signal. Fc receptors specific for IgG, including Fc γ R1 (CD64), Fc γ RII (CD32) and Fc γ RIII (CD16) are present on many cell types, with particularly high expression on myeloid, granulocyte and B cell lineages. Several antibody supply companies provide both human and mouse Fc-blocking reagents that can be used as indicated in the following protocol to minimize non-specific antibody binding.

Reagents and Solutions to Prepare in Advance

Antibody Cocktail

Prepare cocktails of Maxpar Metal Conjugated Antibodies, for both cell surface staining and phospho-proteins, in Cell Staining Buffer. It is recommended to prepare antibody cocktail in a total volume of 50 μ L, so that when added to 50 μ L of cells the total staining volume is 100 μ L. The antibody cocktail can be stored for up to 24 hours before staining.

Intercalation Solution

Prepare 1 mL of cell intercalation solution for each sample by adding Cell-ID Intercalator-Ir into Maxpar Fix and Perm Buffer to a final concentration of 125 nM (a 1000X dilution of the 125 μ M stock solution) and mix by vortexing.

For example, for 10 samples, prepare intercalation solution by adding 10 μL of 125 μM Intercalator-Ir to 10 ml of Fix and Perm Buffer.

[Optional] Serum-Free and Complete Media:

If performing the optional cisplatin viability stain, pre-warm serum-free and complete media at 37°C prior to beginning protocol. Use the same media that is normally used for cell culture.

Protocol

[Optional] Cell-ID Cisplatin Viability Stain (if not performing Cell-ID Cisplatin viability-stain, proceed to step #1 of the protocol below):

Note: The following Cisplatin-staining protocol has minimal impact on the phosphorylation of targets when performed before PMA and Ionomycin cell activation. However, the impact of Cisplatin-staining should be evaluated individually for all phospho-proteins and with different stimulation conditions.

- a. Wash cells to be stimulated with pre-warmed serum-free media and discard supernatant.
 - b. Resuspend cells to 1 X 10⁷/ml in pre-warmed serum-free media and add Cell-ID Cisplatin to final concentration of 5 μM (1000X dilution of 5 mM stock solution).
 - c. Mix well and incubate at 37°C for 5 minutes.
 - d. Quench cisplatin staining by washing with pre-warmed serum-containing complete media using 5-10X the volume of the stained cells, centrifuge and discard supernatant.
 - e. Place cells back in culture conditions for 15 minutes to allow cells to “rest”.
 - f. Proceed with cell activation in step #1 below.
- 1** Prepare cells of interest from cell culture or primary tissue and activate desired signaling pathways by adding stimulus to cells for appropriate length of time.
 - 2** At the end of stimulation, stop the signaling reaction by adding 5X Fix I Buffer to a final concentration of 1X.
 - 3** Mix gently and thoroughly, and incubate for 10 minutes at room temperature.

- 4** Transfer cells to an appropriate tube, and wash with Maxpar Cell Staining Buffer, using 5-10X the volume of the cell suspension; centrifuge and discard supernatant by aspiration.
- 5** Resuspend cells in Maxpar Cell Staining Buffer and aliquot 3 million cells, in a volume of 50 μ L, into 5 mL tubes for each sample to be stained.
- 6** [Optional] Fc-Blocking: add Fc-Receptor Blocking Solution to each tube and incubate for 10 minutes at room temperature; without washing off Fc-Receptor Blocking Solution continue with protocol.
- 7** Add 50 μ L of the surface antibody cocktail to each tube so the total staining volume is 100 μ L.
- 8** Gently vortex samples and incubate for 30 minutes at room temperature.
- 9** Wash by adding 2 mL Maxpar Cell Staining Buffer to each tube, centrifuge and discard supernatant by aspiration.
- 10** Resuspend cells in residual volume by gently vortexing, and place cells on ice for 10 minutes to chill sample.
- 11** Add 1 mL of 4°C methanol to each sample, mix gently, and incubate for 15 minutes on ice.
- 12** Wash cells with 2 mL Maxpar Cell Staining Buffer, centrifuge and discard supernatant by aspiration – repeat for a total 2 washes
- 13** Add phospho-protein antibody cocktail to each tube so the total staining volume is 100 μ L.
- 14** Gently vortex and incubate for 30 minutes at room temperature.
- 15** Following the incubation, wash by adding 2 mL Maxpar Cell Staining Buffer to each tube, centrifuge and discard supernatant by aspiration.
- 16** Repeat for a total of two washes, and resuspend cells in residual volume by gently vortexing after final wash/aspiration.
- 17** Add 1 ml of the intercalation solution to each tube and gently vortex. Incubate for 1 hour at room temperature or leave overnight at 4°C. Note: Cells can be left at 4°C in the intercalation solution up to 48 hours.
- 18** Wash cells by adding 2 ml of Maxpar Cell Staining Buffer, centrifuge and discard supernatant by aspiration.
- 19** Repeat for a total of two washes with Maxpar Cell Staining Buffer.
- 20** Wash cells with 2 ml of Maxpar Water, centrifuge and discard supernatant by aspiration.

21 Leave cells pelleted until ready to run on CyTOF®. Immediately prior to CyTOF data acquisition, adjust cell concentration to $2.5\text{-}5 \times 10^5/\text{ml}$ with Maxpar Water and filter cells into cell strainer cap tubes.

22 Acquire data on CyTOF.

For technical support visit fluidigm.com/support

EMAIL United States and other countries not in Europe or Asia support.northamerica@fluidigm.com | Europe support.europe@fluidigm.com | Asia techsupportasia@fluidigm.com

PHONE United States (toll-free) +1 866 358 4354 | Europe +33 1 60 92 42 40 | Japan +81 3 3662 2150 | China +86 21 3255 8368 | All other countries +1 650 266 6100

For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. Use standard laboratory safety protocols. **Safety data sheet information** fluidigm.com/sds **Patent and license information** fluidigm.com/legalnotices **EU's WEEE directive information** fluidigm.com/compliance | Fluidigm, the Fluidigm logo, Cell-ID, CyTOF, Maxpar are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners. © 2015 Fluidigm Corporation. All rights reserved. 04/2015

Antibody Cocktail Preparation Guide:

The following guide can be used to prepare the Maxpar Metal Conjugated Antibody cocktail in Maxpar Cell Staining Buffer. Prepare the antibody cocktail in a 1.5 ml tube by first adding Cell Staining Buffer and then adding each of the antibodies. Combine 50 µl of the complete antibody cocktail with each sample to be stained.

| (a) Number of Samples | (d) Vol of Antibody (µl) | (b) Number of Antibodies | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--------------------------------|--------------------------------------|--------------------------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|-----|------|------|------|------|------|------|------|------|------|-----|------|------|------|------|
| | | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 | 25 | 26 | 27 | 28 | 29 | 30 | 31 | 32 | 33 | 34 |
| 1 | 1.1 | 53.9 | 52.8 | 51.7 | 50.6 | 49.5 | 48.4 | 47.3 | 46.2 | 45.1 | 44 | 42.9 | 41.8 | 40.7 | 39.6 | 38.5 | 37.4 | 36.3 | 35.2 | 34.1 | 33 | 31.9 | 30.8 | 29.7 | 28.6 | 27.5 | 26.4 | 25.3 | 24.2 | 23.1 | 22 | 20.9 | 19.8 | 18.7 | 17.6 |
| 2 | 2.2 | 108 | 106 | 103 | 101 | 99 | 96.8 | 94.6 | 92.4 | 90.2 | 88 | 85.8 | 83.6 | 81.4 | 79.2 | 77 | 74.8 | 72.6 | 70.4 | 68.2 | 66 | 63.8 | 61.6 | 59.4 | 57.2 | 55 | 52.8 | 50.6 | 48.4 | 46.2 | 44 | 41.8 | 39.6 | 37.4 | 35.2 |
| 3 | 3.3 | 162 | 158 | 155 | 152 | 149 | 145 | 142 | 139 | 135 | 132 | 129 | 125 | 122 | 119 | 116 | 112 | 109 | 106 | 102 | 99 | 95.7 | 92.4 | 89.1 | 85.8 | 82.5 | 79.2 | 75.9 | 72.6 | 69.3 | 66 | 62.7 | 59.4 | 56.1 | 52.8 |
| 4 | 4.4 | 216 | 211 | 207 | 202 | 198 | 194 | 189 | 185 | 180 | 176 | 172 | 167 | 163 | 158 | 154 | 150 | 145 | 141 | 136 | 132 | 128 | 123 | 119 | 114 | 110 | 106 | 101 | 96.8 | 92.4 | 88 | 83.6 | 79.2 | 74.8 | 70.4 |
| 5 | 5.5 | 270 | 264 | 259 | 253 | 248 | 242 | 237 | 231 | 226 | 220 | 215 | 209 | 204 | 198 | 193 | 187 | 182 | 176 | 171 | 165 | 160 | 154 | 149 | 143 | 138 | 132 | 127 | 121 | 116 | 110 | 105 | 99 | 93.5 | 88 |
| 6 | 6.6 | 323 | 317 | 310 | 304 | 297 | 290 | 284 | 277 | 271 | 264 | 257 | 251 | 244 | 238 | 231 | 224 | 218 | 211 | 205 | 198 | 191 | 185 | 178 | 172 | 165 | 158 | 152 | 145 | 139 | 132 | 125 | 119 | 112 | 106 |
| 7 | 7.7 | 377 | 370 | 362 | 354 | 347 | 339 | 331 | 323 | 316 | 308 | 300 | 293 | 285 | 277 | 270 | 262 | 254 | 246 | 239 | 231 | 223 | 216 | 208 | 200 | 193 | 185 | 177 | 169 | 162 | 154 | 146 | 139 | 131 | 123 |
| 8 | 8.8 | 431 | 422 | 414 | 405 | 396 | 387 | 378 | 370 | 361 | 352 | 343 | 334 | 326 | 317 | 308 | 299 | 290 | 282 | 273 | 264 | 255 | 246 | 238 | 229 | 220 | 211 | 202 | 194 | 185 | 176 | 167 | 158 | 150 | 141 |
| 9 | 9.9 | 485 | 475 | 465 | 455 | 446 | 436 | 426 | 416 | 406 | 396 | 386 | 376 | 366 | 356 | 347 | 337 | 327 | 317 | 307 | 297 | 287 | 277 | 267 | 257 | 248 | 238 | 228 | 218 | 208 | 198 | 188 | 178 | 168 | 158 |
| 10 | 11 | 539 | 528 | 517 | 506 | 495 | 484 | 473 | 462 | 451 | 440 | 429 | 418 | 407 | 396 | 385 | 374 | 363 | 352 | 341 | 330 | 319 | 308 | 297 | 286 | 275 | 264 | 253 | 242 | 231 | 220 | 209 | 198 | 187 | 176 |
| 11 | 12.1 | 593 | 581 | 569 | 557 | 545 | 532 | 520 | 508 | 496 | 484 | 472 | 460 | 448 | 436 | 424 | 411 | 399 | 387 | 375 | 363 | 351 | 339 | 327 | 315 | 303 | 290 | 278 | 266 | 254 | 242 | 230 | 218 | 206 | 194 |
| 12 | 13.2 | 647 | 634 | 620 | 607 | 594 | 581 | 568 | 554 | 541 | 528 | 515 | 502 | 488 | 475 | 462 | 449 | 436 | 422 | 409 | 396 | 383 | 370 | 356 | 343 | 330 | 317 | 304 | 290 | 277 | 264 | 251 | 238 | 224 | 211 |
| 13 | 14.3 | 701 | 686 | 672 | 658 | 644 | 629 | 615 | 601 | 586 | 572 | 558 | 543 | 529 | 515 | 501 | 486 | 472 | 458 | 443 | 429 | 415 | 400 | 386 | 372 | 358 | 343 | 329 | 315 | 300 | 286 | 272 | 257 | 243 | 229 |
| 14 | 15.4 | 755 | 739 | 724 | 708 | 693 | 678 | 662 | 647 | 631 | 616 | 601 | 585 | 570 | 554 | 539 | 524 | 508 | 493 | 477 | 462 | 447 | 431 | 416 | 400 | 385 | 370 | 354 | 339 | 323 | 308 | 293 | 277 | 262 | 246 |
| 15 | 16.5 | 809 | 792 | 776 | 759 | 743 | 726 | 710 | 693 | 677 | 660 | 644 | 627 | 611 | 594 | 578 | 561 | 545 | 528 | 512 | 495 | 479 | 462 | 446 | 429 | 413 | 396 | 380 | 363 | 347 | 330 | 314 | 297 | 281 | 264 |
| 16 | 17.6 | 862 | 845 | 827 | 810 | 792 | 774 | 757 | 739 | 722 | 704 | 686 | 669 | 651 | 634 | 616 | 598 | 581 | 563 | 546 | 528 | 510 | 493 | 475 | 458 | 440 | 422 | 405 | 387 | 370 | 352 | 334 | 317 | 299 | 282 |
| 17 | 18.7 | 916 | 898 | 879 | 860 | 842 | 823 | 804 | 785 | 767 | 748 | 729 | 711 | 692 | 673 | 655 | 636 | 617 | 598 | 580 | 561 | 542 | 524 | 505 | 486 | 468 | 449 | 430 | 411 | 393 | 374 | 355 | 337 | 318 | 299 |
| 18 | 19.8 | 970 | 950 | 931 | 911 | 891 | 871 | 851 | 832 | 812 | 792 | 772 | 752 | 733 | 713 | 693 | 673 | 653 | 634 | 614 | 594 | 574 | 554 | 535 | 515 | 495 | 475 | 455 | 436 | 416 | 396 | 376 | 356 | 337 | 317 |
| 19 | 20.9 | 1024 | 1003 | 982 | 961 | 941 | 920 | 899 | 878 | 857 | 836 | 815 | 794 | 773 | 752 | 732 | 711 | 690 | 669 | 648 | 627 | 606 | 585 | 564 | 543 | 523 | 502 | 481 | 460 | 439 | 418 | 397 | 376 | 355 | 334 |
| 20 | 22 | 1078 | 1056 | 1034 | 1012 | 990 | 968 | 946 | 924 | 902 | 880 | 858 | 836 | 814 | 792 | 770 | 748 | 726 | 704 | 682 | 660 | 638 | 616 | 594 | 572 | 550 | 528 | 506 | 484 | 462 | 440 | 418 | 396 | 374 | 352 |
| 21 | 23.1 | 1132 | 1109 | 1086 | 1063 | 1040 | 1016 | 993 | 970 | 947 | 924 | 901 | 878 | 855 | 832 | 809 | 785 | 762 | 739 | 716 | 693 | 670 | 647 | 624 | 601 | 578 | 554 | 531 | 508 | 485 | 462 | 439 | 416 | 393 | 370 |
| 22 | 24.2 | 1186 | 1162 | 1137 | 1113 | 1089 | 1065 | 1041 | 1016 | 992 | 968 | 944 | 920 | 895 | 871 | 847 | 823 | 799 | 774 | 750 | 726 | 702 | 678 | 653 | 629 | 605 | 581 | 557 | 532 | 508 | 484 | 460 | 436 | 411 | 387 |
| 23 | 25.3 | 1240 | 1214 | 1189 | 1164 | 1139 | 1113 | 1088 | 1063 | 1037 | 1012 | 987 | 961 | 936 | 911 | 886 | 860 | 835 | 810 | 784 | 759 | 734 | 708 | 683 | 658 | 633 | 607 | 582 | 557 | 531 | 506 | 481 | 455 | 430 | 405 |
| 24 | 26.4 | 1294 | 1267 | 1241 | 1214 | 1188 | 1162 | 1135 | 1109 | 1082 | 1056 | 1030 | 1003 | 977 | 950 | 924 | 898 | 871 | 845 | 818 | 792 | 766 | 739 | 713 | 686 | 660 | 634 | 607 | 581 | 554 | 528 | 502 | 475 | 449 | 422 |
| 25 | 27.5 | 1348 | 1320 | 1293 | 1265 | 1238 | 1210 | 1183 | 1155 | 1128 | 1100 | 1073 | 1045 | 1018 | 990 | 963 | 935 | 908 | 880 | 853 | 825 | 798 | 770 | 743 | 715 | 688 | 660 | 633 | 605 | 578 | 550 | 523 | 495 | 468 | 440 |

TO USE THE TABLE: Locate the row matching the number of samples to be processed (a) and the column for the number of antibodies used to stain the sample (b). Use the table to determine the Total Volume of Cell Staining Buffer needed (c). Add this volume of Cell Staining Buffer to your mastermix tube. Again locate the row matching the number of samples to be processed (a) and in the adjacent column determine the volume per antibody (d). Add the indicated volume of each antibody solution to the mastermix tube.