



BD FACSymphony™ A1 Cell Analyzer

Premium performance in a benchtop footprint





BD FACSymphony™ A1 Cell Analyzer features:


- Premium high-end BD FACSymphony™ instrument technology scaled to fit on your benchtop
- Flexibility to meet a broad spectrum of research needs from **small particle research** to **16-color immunophenotyping**
- Industry standard **BD FACSDiva™ Software** for streamlined workflow from system setup to data acquisition and analysis

BD® Small Particle Detector

Optional detector for small particle side scatter (SP SSC) for resolving particles as small as 90 nm

Familiar Workflow

Industry standard BD FACSDiva™ and FlowJo™ Software enable easy system QC, data acquisition and analysis



Trusted Partner
Supported by over 45 years of flow cytometry expertise

BD FACSymphony™ Instrument Technology

Reduce background noise and increase sensitivity with low-noise electronics, high-powered lasers and tight beam spot, shared across the premium BD FACSymphony™ platforms. The BD FACSymphony™ A1 Cell Analyzer is compatible with BD Horizon™ Dyes and supports up to 16 colors or 19 parameters simultaneously.





Premium BD FACSymphony™ instrument technology

delivered in a compact size

Up to

16 fluorochromes and 19 parameters

to conduct deep and broad phenotyping



Enhance detection sensitivity

with four high-powered 100 mW lasers: Violet (405 nm), Blue (488 nm), Yellow-Green (561 nm) and Red (637 nm)

BD® Small Particle Detector Option for analysis of small particles

such as extracellular vesicles including exosomes



Identify and analyze rare cell types and events

with our redesigned optics including small beam spots combined with low-noise electronics



Gain rich scientific insights

by leveraging BD Horizon Brilliant™ Reagents



Enable easy system QC

using industry-standard BD FACSDiva™ Software and BD® CS&T Beads



Utilizes FlowJo™ Software, the leading bioinformatics platform*

for flow cytometry analysis



Ideal for labs with limited space

Small footprint (58 x 61 x 59 cm)

Automated sample processing in high-throughput mode

using the BD® High-Throughput Sampler Option



*In 2021, FlowJo™ Software was cited in leading immunology peer-reviewed journals 80% of the time a flow cytometry analysis software package was cited.

Able to detect up to 16 colors and resolve rare cell subsets

Table 1. Instrument configuration and reagents in the cytotoxic immune cells panel

Laser	Filter	Fluorochrome	Specificity
Violet 405 nm	450/50	BV421	Perforin
	525/50	BV480	CD159a (NKG2A)
	610/20	BV605	CD19
			CD14
			CD123
			CD141
		FV5575V	-
	670/30	BV650	CD3
710/50	BV711	CD314 (NKG2D)	
780/60	BV786	HLA-DR	
Blue 488 nm	530/30	FITC	CD57
	710/50	PerCP-Cy5.5	CD8
Yellow-Green 561 nm	586/15	PE	CD158 (KIRs)
	610/20	PE-CF594	CD56
	670/30	PE-Cy5	CD95 (Fas)
	710/50	PE-Cy5.5	CD127 (IL7R- α)
	780/60	PE-Cy7	CD38
Red 637 nm	670/30	AF647	Granzyme K
	710/50	R718	Granzyme B
	780/60	APC-H7	CD16 (Fc γ RIII)

BV, BD Horizon Brilliant™ Violet; FVS, BD Horizon™ Fixable Viability Stain; BD Horizon™ Red, Alexa Fluor™

Figure 1




Figure 1. Identification of cytotoxic immune cell populations in healthy human peripheral blood

Within live and lineage negative cells, analysis of CD56 versus CD3 revealed various cell populations that were color coded as cytokine-producing NK cells (pink), cytotoxic NK cells (blue), CD56⁺ T cells containing NKT cells (brown), CD56⁺ T cells containing $\gamma\delta$ T cells (green) and cytotoxic CD8⁺ T cells (purple).

To learn more, download the panel sheet *Characterization of Cytotoxic Immune Cells in Human Peripheral Whole Blood* from bdbiosciences.com

Figure 2A




Figure 2B






Figure 2. Phenotyping of circulating cytotoxic cells using a 16-color panel


The plots represent the analysis of cytolytic proteins in combination with various cell differentiation markers, enabling a deeper characterization of the cell populations gated in Figure 1. A. Overlay of NK cell subsets. B. Identification of activated CD8 T cells based on the expression of CD38 and HLA-DR. The HLA-DR FMO staining helped to determine the gating boundaries for proper detection of the double positive cells.

Independent detection of large (SSC) and small (SP SSC) particles

The BD FACSymphony™ A1 Cell Analyzer with optional BD® Small Particle Detector is able to resolve scatter of small particles such as extracellular vesicles, viral particle, exosomes and more.



Seamless small particle detection workflow





*The BD® Small Particle Detector can resolve particles as small as 90 nm.


†In 2021, FlowJo™ Software was cited in leading immunology peer-reviewed journals 80% of the time a flow cytometry analysis software package was cited.

Detection of extracellular vesicles using BD FACSymphony™ A1 Cell Analyzer



Resolution of 90-nm polystyrene particles with the BD® Small Particle Detector option



Side scatter sensitivity of BD® Small Particle Detector



Characterization of extracellular vesicles from human MCF7 cell line*



* This work was performed in collaboration with Wauben Lab (Utrecht University, The Netherlands) and was supported by the TRAIN-EV Marie Skłodowska-Curie Action-Innovative Training Network, <http://train-ev.eu> grant agreement No 722148


BD FACSymphony™ Systems



	BD FACSymphony™ A1	BD FACSymphony™ A3	BD FACSymphony™ A5	BD FACSymphony™ A5 SE
Number of lasers	4	5	5–9	5
Fluorescent detectors	16	Up to 28	Up to 48	Up to 48
Instrument type	Analyzer	Analyzer	Analyzer	Analyzer
Software	BD FACSDiva™	BD FACSDiva™	BD FACSDiva™	BD FACSDiva™
Footprint	58 x 61 cm	83.8 x 76.2 cm	101.6 x 78.7 cm	101.6 x 78.7 cm
Small particle detector	Yes	Custom	Custom	Custom

Backed and supported by BD

We're committed to partnering with you to provide the mission-critical tools and support you need to advance your research.



BD Life Sciences – Biosciences Regional Offices

bdbiosciences.com/contact

Australia

Toll Free 1800.656.100
Tel 61.2.8875.7000
Fax 61.2.8875.7200

Canada

Tel 866.979.9408
Fax 888.229.9918

China

Tel 86.21.3210.4610
Fax 86.21.5292.5191

Europe

Tel 32.2.400.98.95
Fax 32.2.401.70.94

India

Tel 91.124.2383566
Fax 91.124.2383224/25/26

Japan

Nippon Becton Dickinson
Toll Free 0120.8555.90
Fax 81.24.593.3281

Latin America/Caribbean

Toll Free 0800.771.71.57
Tel 55.11.5185.9688

New Zealand

Toll Free 0800.572.468
Tel 64.9.574.2468
Fax 64.9.574.2469

Singapore

Tel 65.6690.8691
Fax 65.6860.1593

United States

US Orders 855.236.2772
Technical Service
877.232.8995
Fax 800.325.9637

Office locations are available on our websites.

Class 1 Laser Product.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

BD Life Sciences, San Jose, CA 95131, USA

bdbiosciences.com

BD, the BD Logo, BD FACSDiva, BD FACSymphony, BD Horizon Brilliant Violet, BD Horizon Fixable Viability Stain and FlowJo are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners.
© 2022 BD. All rights reserved. BD-30354 (v3.0) 0722

Alexa Fluor (AF) is a trademark of Thermo Fisher Scientific. Cy is a trademark of Global Life Sciences Germany GmbH or an affiliate doing business as Cytiva.

