

# Cell-ID™ Intercalator-Ir 500 µM

Catalog #201192B (500 µL)

**WARNING! CHEMICAL HAZARD.** Before handling any chemicals, refer to the Safety Data Sheet (SDS) provided by the manufacturer, and observe all relevant precautions.

**NOTICE: HIGH CONCENTRATION.** Cell-ID Intercalator-Ir 500 µM is a highly concentrated metal intercalator solution and must be diluted in accordance with this protocol to avoid early failure of the detector.

## Description

Cell-ID Intercalator-Ir is a cationic nucleic acid intercalator that contains natural abundance Iridium (<sup>191</sup>Ir and <sup>193</sup>Ir) and is used for identifying nucleated cells in CyTOF® system analysis. When cells are stained with Intercalator-Ir, it will bind to cellular nucleic acid, and detection of both stable isotopes will enable identification of nucleated cells. It is a live cell membrane-impermeable dye and therefore requires cells to be fixed and/or permeabilized before staining.

**Note:** While dilutions of the 500 µM stock solution are suggested in the protocols below, the concentration can be titrated for individual cell types and experiments for optimal Cell-ID Intercalator staining. It is suggested not to exceed 1 µM intercalator concentration in the staining solution.

## Storage

Upon arrival, aliquot and store at -20 °C. Aliquots stored at 4 °C are stable for up to three months.

## Staining Protocol

- 1 After cell staining is complete, prepare 1 mL of cell intercalation solution for each sample by diluting Cell-ID Intercalator-Ir 1:4000 into Maxpar® fix and perm buffer (Fluidigm Cat. 201067) and mix by vortexing.
- 2 Add 1 mL of the intercalation solution prepared in step 1 to each tube and gently vortex. Incubate for one hour at room temperature or leave overnight at 4 °C.

**Note:** Cells can be left at 4 °C in the intercalation solution up to 48 hours.

- 3 Wash cells by adding 2 mL of Maxpar® cell staining buffer (Fluidigm Cat. 201068), centrifuge and discard supernatant by aspiration.
- 4 Repeat for a total of two washes with Maxpar cell staining buffer.
- 5 Wash cells with 2 mL of Maxpar® water (Fluidigm Cat. 201069), centrifuge and discard supernatant by aspiration.
- 6 Leave cells pelleted until ready to run on the CyTOF system. Immediately prior to data acquisition, adjust cell concentration to 2.5–5 x 10<sup>5</sup>/mL with Maxpar water and filter cells into cell strainer cap tubes.
- 7 Acquire data on the CyTOF system.

**For technical support visit [fluidigm.com/support](http://fluidigm.com/support)**

**North America** +1 650 266 6100 | Toll-free: +1 866 358 4354 in the US | [support.northamerica@fluidigm.com](mailto:support.northamerica@fluidigm.com) **Europe** +33 1 60 92 42 40 | [support.europe@fluidigm.com](mailto:support.europe@fluidigm.com)  
**China (excluding Hong Kong)** +86 21 3255 8368 | [techsupportchina@fluidigm.com](mailto:techsupportchina@fluidigm.com) **Japan** +81 3 3662 2150 | [techsupportjapan@fluidigm.com](mailto:techsupportjapan@fluidigm.com)  
**All other Asian countries** +1 650 266 6100 | [techsupportasia@fluidigm.com](mailto:techsupportasia@fluidigm.com) **Central and South America** +1 650 266 6100 | [techsupportlatam@fluidigm.com](mailto:techsupportlatam@fluidigm.com)

**For Research Use Only. Not for use in diagnostic procedures.**

Information in this publication is subject to change without notice. **Safety data sheet information** [fluidigm.com/sds](http://fluidigm.com/sds) **Patent and license information** [fluidigm.com/legalnotices](http://fluidigm.com/legalnotices) | Fluidigm, the Fluidigm logo, Cell-ID, CyTOF, and Maxpar are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2015 Fluidigm Corporation. All rights reserved. 10/2015